

ECE Group 0/1
0 - 18 kg
0 - 39.7 lb

CO-PILOT

The lightest rehab reboarder



THOMASHILFEN Co-Pilot

only 5 kg / 11 lb




LIGHTWEIGHT

Back-friendly everyday handling for parents and carers – no heavy carrying when seat is used in several vehicles.

UPHOLSTERY

Cosy upholstery, the breathable cover can be completely removed and washed at 40° C.



SHOULDER HEIGHT

The padded belt grows with the child: 3-fold height adjustment possible



REHAB EQUIPMENT

- including:
- Abduction block
 - Hip supports
 - Head pad

BELT LENGTH

Quick and easy adjustment of the belt length via central adjuster at the front of the seat



“ With the lightest rear facing rehab car seat (up to 10 kg / 22 lb body weight), everyday handling becomes child’s play. And the delicate child’s head is safely supported and protected.

When the children grow bigger (from 9 kg / 19.8 lb), the THOMASHILFEN Co-Pilot is simply turned around and the direction of travel is changed. In the different seat positions the child can be active, rest or sleep – whatever the situation requires. ”

AT A GLANCE

The current seating position can be quickly read off the markings:

- ➔ blue: rear-facing
- ➔ red: 3 positions forward-facing

CHANGING POSITION

With one quick flick of the hand, the forward facing car seat can be adjusted to 3 comfortable sitting and resting positions



Thomashilfen

In everyday use



Rear facing (up to 10 kg / 22 lb)

as a safe reboarder – the child’s head is perfectly protected in the event of a collision



Forward facing

(from 9 to 18 kg / 19.8 to 39.7 lb)

car seat can be used in 3 comfortable sitting and resting positions



Uncomplicated installation of the reboarder in the car – without annoying, long assembly



YouTube

Take a look at the video:
How to install the car seat
THOMASHILFEN Co-Pilot in a car

THOMASHILFEN Co-Pilot



All details online:
www.thomashilfen.com/co-pilot



Features



Removable abduction block for a good sitting position, prevents the legs from crossing over



Thick high hip supports with adjustable padding thickness provide extra guidance and stability for smaller children



Head pad can be precisely attached for selective support of the sensitive child's head.



Padded 5-point safety belt with crotch pad and central adjuster for quick adjustment of the belt length

Accessory

5 upholstery designs available



Lateral trunk pads can be precisely attached to the side of the seat for better upper body guidance



GREY

RED

BLUE

GREEN

MAGENTA

Soft, breathable cover for pleasant sitting comfort without annoying sweating. The covers are removable and washable at 40° C. Do not tumble dry.

You will find complete overview of all accessories on the order form or our website:
www.thomashilfen.com

Dimensions & weights	THOMASHILFEN Co-Pilot
Seat depth	29.5 cm / 11.6"
Seat width front	23.5 cm / 9.3" (20,5 cm / 8.1" with hip supports)
Seat width back	20 cm / 7.9" (14.5 cm / 5.7" with hip supports)
Back height (seat)	51 cm / 20.1"
Shoulder height (belt guide)	19 / 25 / 31 cm // 7.5 / 9.8 / 12.2"
Seat tilt-in-space (rear-facing)	45°
(forward-facing)	0° / 10° / 20°
Overall dimensions (w x h x d)	42 x 62,5 x 43 cm / 16.5 x 24.6 x 16.9"
Weight of seat	ca. 5 kg / 11 lb
max. user weight	rear-facing: 0 to 10 kg / 0 to 22 lb forward-facing: 9 to 18 kg / 19.8 to 39.7 lb

Basic configuration:

Seat shell (rearward facing: 1 seating position / forward facing: 3 seating positions); seat cushion; head pad; crotch pad; 5-point harness with central adjuster; abduction block; hip supports (coloured, matching the basic model); removable covers, machine washable up to 40° C; do not tumble dry

Co-Pilot

Rückhalteeinrichtung für spezielle Anforderungen
 Restraint system for special needs

S
 0 – 18 kg
 0 – 39.7 lb



Thomas Hilfen
 GmbH & Co. Medico KG
 27432 Bremervörde
 Germany



Thomas Hilfen für Körperbehinderte GmbH & Co. Medico KG · Walkmühlenstraße 1 · D - 27432 Bremervörde
 Phone: +49 (0)4761 886-63 or -68 · Fax: +49 (0)4761 886-19 · www.thomashilfen.com