



Tendercare Ltd
Unit 10, Minster Court,
Courtwick Lane,
Littlehampton
BN17 7RN

Tel 01903 726161
Email info@tendercareltd.com
Web www.tendercareltd.com

Child's Name.....

Contact Name.....

Address.....

.....

.....

Tel No.....

Email

Hercules Small R129 Car Seat

TH470 Hercules Small R129 (user weight 18-50kg)

Select the side of pre-assembled belt guide (UK car)

Behind driver

Behind Passenger

Grey/Blue 

Grey/Red 

✓ TH479 ISOFIX/Seatfix Connection (Hercules Small R129)

Please note, the car's own seat belt must be used in front of the occupant at all time. Isofix connection is required to fit the car seat /swivel/tilting base.

Hercules Small R129 Car Seat

ECE R129/03



Technical Information

Seat Depth	310-420mm
Seat Width	290-380mm
Back Height	570-800mm
Shoulder Width	420mm
Shoulder Height (belt guide)	330-600mm
Lower Leg Length (with footrest)	160-420mm
Max. External Dims. (w x h x d)	550 x 850 x 340mm
Min. Folded Dims. (w x h x d)	550 x 230 x 600mm
Weight of Seat	11.6kg
User Weight Capacity	18-50kg

Hercules Small Car Seat includes:

Seat and back with integrated seat depth and back height adjustments, continuous back angle adjustment with quick release lever, breathable upholstery with flame-retardant cover fabric, lateral headrest pads, back cushion, hip & leg guides, seat cushion, crotch pad, 5-point belt including strap soft pads, 1 x belt guide for 3-point vehicle belt, side impact protectors, prepared mounting options for all accessories, removable covers - machine washable up to 40°C.

OPTIONS FOR FURTHER STABILITY

TH480 Front Stabilizing Stand - recommended for stability if footrest not being fitted.
(Stand and footrest cannot be used in combination)



TH513 Cushion for swivel/tilting base (recommended for vehicle seats with seat recess to compensate for the contact surface of the function plate and vehicle seat)



ACCESSORIES

TH4751 Swivel Base

(For better support we recommend to use the swivel base with stand and/or cushion)

Please select the side of swivel function belt:

BEHIND DRIVER

BEHIND PASSENGER



TH4761 Tilting Base 10° (0°/2.5°/5°/7.5°/10°) PRE-MOUNTED - Must not be adjusted when client is in seat.

(For better support we recommend to use the tilting function with stand and/or cushion)



TH477 Tilting Base Extension 12.5° to 20° (12.5°/15°/17.5°/20°)
- only in combination with Tilting function 10°




OPTIONAL ACCESSORIES

TH494 Head additional cushion



TH483 Shoulder pads (Pair) 4cm each



TH486 Lateral Trunk Support - Firm



TH487 Lateral Trunk Support - Swing away




Only one set of Lateral Trunk Supports can be selected.


TH498 Crotch Pad with magnetic lock
(to help prevent unlocking of harness buckle)



TH488 Pommel/Abduction block (height adjustable)



TH490 Additional seat belt guide, top



TH423 Chest Clip - Fits to 5 point harness



TH489 Armrests



TH491 Table (armrests are required for fitting)



TH481 Footrest, foldable (height and angle adjustable.)



TH8141 Car seat protection mat including back - requires vehicle's headrest poles for fixing.



TH497 Incontinence Seat Cover
(fits on top of seat base cushion)

It is possible to insert other commercially available absorbent pads into the incontinence seat pad to increase the absorption capacity if required.

ADDITIONAL INFORMATION

Height of childcm

Weight of childkg

Any Additional Information:

.....

.....

When seated:
Base of bottom to top of head.....cm
Base of bottom to top of shoulder.....cm