

 Tendercare Tendercare Ltd, Unit 10, Minster Court, Littlehampton, West Sussex BN17 7RN Tel 01903 726161 Email: info@tendercareltd.com www.tendercareltd.com	Octopus Positioning Car Harness The Octopus Positioning Harness is not classified as a safety harness, but as a Medical Device for Positioning.  <table border="1" style="width: 100%; border-collapse: collapse; text-align: center;"> <thead> <tr> <th style="padding: 5px;">Technical Information</th> <th style="padding: 5px;">SMALL</th> <th style="padding: 5px;">MEDIUM</th> </tr> </thead> <tbody> <tr> <td style="padding: 5px;">7 Point Seat Vest W x H</td> <td style="padding: 5px;">33x41 cm 13"x16.1"</td> <td style="padding: 5px;">42x42 cm 16.5"x16.5"</td> </tr> <tr> <td style="padding: 5px;">Harness Base Back Part W x H</td> <td style="padding: 5px;">37x82 cm 14.6"x32.3"</td> <td style="padding: 5px;">37x82 cm 14.6"x32.3"</td> </tr> <tr> <td style="padding: 5px;">Body height user</td> <td style="padding: 5px;">130-160cm 51.2"-63"</td> <td style="padding: 5px;">150-190cm 59.1"-74.8"</td> </tr> <tr> <td style="padding: 5px;">Body weight user</td> <td style="padding: 5px;">30-50 kg 66.1-110.2 lbs</td> <td style="padding: 5px;">50-80 kg 110.2-176.3 lbs</td> </tr> </tbody> </table>	Technical Information	SMALL	MEDIUM	7 Point Seat Vest W x H	33x41 cm 13"x16.1"	42x42 cm 16.5"x16.5"	Harness Base Back Part W x H	37x82 cm 14.6"x32.3"	37x82 cm 14.6"x32.3"	Body height user	130-160cm 51.2"-63"	150-190cm 59.1"-74.8"	Body weight user	30-50 kg 66.1-110.2 lbs	50-80 kg 110.2-176.3 lbs
Technical Information	SMALL	MEDIUM														
7 Point Seat Vest W x H	33x41 cm 13"x16.1"	42x42 cm 16.5"x16.5"														
Harness Base Back Part W x H	37x82 cm 14.6"x32.3"	37x82 cm 14.6"x32.3"														
Body height user	130-160cm 51.2"-63"	150-190cm 59.1"-74.8"														
Body weight user	30-50 kg 66.1-110.2 lbs	50-80 kg 110.2-176.3 lbs														
Client Name..... Contact Name..... Address..... Tel No..... Email	Octopus Positioning Car Harness <i>Weight and heights are for guidance only. Correct size to be chosen to ensure comfort and fitting for the user.</i> TH693 Small - (user weight 30-50 kg) - (user height 130-160 cm) TH694 Medium - (user weight 50-80 kg) - (user height 150-190 cm) For heights under 135 cm or up to the age of 12, the Octopus Positioning Harness must be used in conjunction with a booster seat. For ease of fitting place the booster seat (NOT High Back) on the car seat after the Octopus Harness has been mounted. The vehicle 3-point seat belt must always be worn over the Octopus positioning harness, regardless of height.															
OPTIONAL ACCESSORIES																
TH691 Belt extension set - If the vehicle seat does not have an opening between the seat and backrest you will need the belt extension accessory to use the product safely.  ONLY required when using one harness in two vehicles - base back parts remain completely mounted in the vehicle. TH695 Basic back part set - base back part, chest strap, shoulder straps, lap belts, crotch strap and breathable padding. 	TH80627 Headrest - Provides lateral support for the head, with height-adjustable support pads.  TH698 Seat Cushion TH699 Incontinence Cover - for Octopus seat cushion 															
ADDITIONAL INFORMATION																
Height of childcm Weight of childkg Any Additional Information:	When seated: Base of bottom to top of head.....cm Base of bottom to top of shoulder.....cm															